

Lab Report 1 - BIOINSTRUMENTATION

Group 11:

Yi Shao Liao

Xuchuan Wang

Rentao Song


Experiment description:

To test the sensing performance of L-glutamate sensor, we add different solution with different concentration of L-glutamate and also different composition at different time. In the end we verify the performance by sensor's sensitivity and limit of detection (LOD).


Experiment setup:


We record the current drop of L-glutamate sensing electrode and Self-referencing for L-glutamate electrode with different solution adding at time showing below.

Time	Added Solution
50s	20 μ M Ascorbic Acid
100s	10 μ M L-glutamate
150s	30 μ M L-glutamate
200s	50 μ M L-glutamate
260s	70 μ M L-glutamate
310s	90 μ M L-glutamate
350s	2 μ M Dopamine
400s	Hydrogen Peroxide



Results:





$$\text{Sensitivity} = 14.7 \left(\frac{\text{nA}}{\text{mm}^2 \times \mu\text{M}} \right)$$

$$\text{Limit of Detection (LOD)} = \frac{(3 \times \text{standard deviation of the baseline})}{(\text{least squares slope})}$$

$$\text{LOD} = 11.26 (\mu\text{M})$$